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LITERATUR	
IFCC	„Approved IFCC Recommendations on Whole Blood Sampling, Transport and Storage for Simultaneous etermination of pH, Blood Gases and Electrolytes“ Eur J Clin Chem Clin Biochem. 1995 Apr;33(4):247-53.
AARC	Clinical Practice Guideline “Blood Gas Analysis and Hemoximetry: 2013”, published on July 30, 2013 as DOI: 10.4187/respcare.02786
AARC	Clinical Practice Guideline “Capillary Blood Gas Sampling for Neonatal & Pediatric Patients” - Respiratory Care (Respir Care 2001;46(5):506–513)
AARC	Clinical Practice Guideline “Sampling for Arterial Blood Gas Analysis” - Respiratory Care (Respir Care 1992;37:913-917)

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EN	Plastic Capillaries for Blood Gas
PRODUCT DESCRIPTION	
Intended use:	In vitro diagnostic medical device for capillary blood sampling with or without an anticoagulant, sample container for the collection and dispensing of human blood samples for subsequent pH/ blood gas, oximetry, electrolyte and metabolite analyses
Filling volume:	See packaging
Anticoagulant:	Amount in I.U.: See packaging
Storage conditions:	2 °C - 30 °C / store in a dry place / protect from sunlight

INSTRUCTIONS FOR USE AND RESTRICTIONS

- The plastic capillaries for blood gas may only be used by trained medical personnel.
- Adhere to the applicable hygiene regulations when handling biological sample material. It may contain dangerous pathogens. Protective gloves, safety glasses and suitable protective clothing must be worn. Risk of infection!
- The analysis should be carried out immediately after collecting the sample. If the sample is stored in the plastic capillaries for blood gas before measurement, sedimentation may occur and this may affect the measured values. This is especially the case for samples with high haematocrit levels, extreme measurement values, elevated auto-metabolism, or increased temperature!
- This product is not suitable for invasive applications!

Please refer to the relevant literature.

The plastic capillaries for blood gas are are a Single use product and not suitable for reuse and must be properly disposed of after use.

BLLOOD SAMPLE COLLECTION: Without mixing sticks, Six-edge closing caps, and mixing magnet:

- Decisions relating to the puncture site and collection (including skin disinfection) should be made by the physician carrying out the treatment.
- Stimulate blood flow in the puncture site area, for example, using a heat pack.
- Perforate the skin so that drops of blood form quickly. Discard the first drop of blood. Fill the plastic capillary for blood gas with blood from the middle of the drop of blood without the formation of air bubbles.
- Aspirate the sample immediately in the analysis device or blood gas analysing system.**

BLOOD SAMPLE COLLECTION With mixing sticks, Six-edge closing caps, and mixing magnets.

- Decisions relating to the puncture site and collection (including skin disinfection) should be made by the physician carrying out the treatment.
- Stimulate blood flow in the puncture site area, for example, using a heat pack.
- Loosely place a six-edge closing cap on the end of the plastic capillary for blood gas.
- Insert a mixing stick into the plastic capillary for blood gas and allow it to slide until it reaches the loosely placed six-edge closing cap.
- Perforate the skin so that drops of blood form quickly. Discard the first few drops of blood. Fill the plastic capillary for blood gas with blood from the middle of the drop of blood without the formation of air bubbles.
- Firmly attach the six-edge closing cap. Then firmly attach a second six-edge closing cap at the opposing end of the capillary.
- Using the mixing magnet, move the mixing stick to and fro along the entire length of the capillary 20 times to mix the blood with the anticoagulant.
- Analyse the sample within 10 minutes. If necessary, store the sample horizontally at ~ 4 °C for no longer than 20 minutes.
- Mix the sample before carrying out the analysis. Do not remove the mixing stick and slide it to the end of the capillary through which no blood will be aspirated.
- Remove both six-edge caps.
- Immediately aspirate the sample into the blood gas/ISE analysis device.**

You can find further information in the blood gas/ISE analysis device manual. For optimum analysis results, use suitable accessories, such as mixing sticks, Six-edge caps, capillary luer adaptor with blood clot catcher and mixing magnet, from our product range.

DISPOSAL

Dispose of the plastic capillaries for blood gas in accordance with the applicable local or laboratory regulations (biological contaminated - hazardous waste)

LITERATURE	
IFCC	“Approved IFCC Recommendations on Whole Blood Sampling, Transport and Storage for Simultaneous Determination of pH, Blood Gases and Electrolytes“ Eur J Clin Chem Clin Biochem. 1995 Apr;33(4):247-53.
AARC	Clinical Practice Guideline “Blood Gas Analysis and Hemoximetry: 2013”, published on July 30, 2013 as DOI: 10.4187/respcare.02786
AARC	Clinical Practice Guideline “Capillary Blood Gas Sampling for Neonatal & Pediatric Patients” - Respiratory Care (Respir Care 2001;46(5):506–513)
AARC	Clinical Practice Guideline “Sampling for Arterial Blood Gas Analysis” - Respiratory Care (Respir Care 1992;37:913-917)

	This product fulfils the requirements of the German Medical Devices Law and European Directive 98/79/EC for in vitro diagnostic medical devices
	In vitro diagnostic medical device
	For single blood collection only. Do not reuse or refurbish
	Store at 2 °C - 30 °C
	Store in a cool dry place
	Keep cool. Protect from sunlight
	Manufacturer: SC – Sanguis Counting Kontrollblutherstellungs- und Vertriebs GmbH Im Löchelchen 11 – 51588 Nümbrecht – Germany www.sanguis-counting.de/ifu/pcbfg sales@sanguis-counting.de +49 2293 9129 0

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
Referencias	
IFCC	«Approved IFCC Recommendations on Whole Blood Sampling, Transport and Storage for Simultaneous Determination of pH, Blood Gases and Electrolytes» Eur J Clin Chem Clin Biochem. 1995 Apr;33(4):247-53.
AARC	Clinical Pracice Guideline «Blood Gas Analysis and Hemoximetry: 2013», published on July 30, 2013 as DOI: 10.4187/respcare.02786
AARC	Clinical Practice Guideline «Capillary Blood Gas Sampling for Neonatal & Pediatric Patients” - Respiratory Care (Respir Care 2001;46(5):506–513)
AARC	Clinical Practice Guideline «Sampling for Arterial Blood Gas Analysis” - Respiratory Care (Respir Care 1992;37:913-917)

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
Bibliographie	
IFCC	“Approved IFCC Recommendations on Whole Blood Sampling, Transport and Storage for Simultaneous Determination of pH, Blood Gases and Electrolytes” Eur J Clin Chem Clin Biochem. 1995 Apr;33(4):247-53.
AARC	Clinical Practice Guideline “Blood Gas Analysis and Hemoximetry: 2013”, published on July 30, 2013 as DOI: 10.4187/respcare.02786
AARC	Clinical Practice Guideline “Capillary Blood Gas Sampling for Neonatal & Pediatric Patients” - Respiratory Care (Respir Care 2001;46(5):506–513)
AARC	Clinical Practice Guideline “Sampling for Arterial Blood Gas Analysis” - Respiratory Care (Respir Care 1992;37:913-917)

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
Letteratura	
IFCC	“Approved IFCC Recommendations on Whole Blood Sampling, Transport and Storage for Simultaneous Determination of pH, Blood Gases and Electrolytes” Eur J Clin Chem Clin Biochem. 1995 Apr;33(4):247-53.
AARC	Clinical Practice Guideline “Blood Gas Analysis and Hemoximetry: 2013”, published on July 30, 2013 as DOI: 10.4187/respcare.02786
AARC	Clinical Practice Guideline “Capillary Blood Gas Sampling for Neonatal & Pediatric Patients” - Respiratory Care (Respir Care 2001;46(5):506–513)
AARC	Clinical Practice Guideline “Sampling for Arterial Blood Gas Analysis” - Respiratory Care (Respir Care 1992;37:913-917)

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Literatura	
IFCC	Approved IFCC Recommendations on Whole Blood Sampling, Transport and Storage for Simultaneous Determination of pH, Blood Gases and Electrolytes" Eur J Clin Chem Clin Biochem. 1995 Apr;33(4):247-53.
AARC	Clinical Practice Guideline „Blood Gas Analysis and Hemoximetry: 2013", published on July 30, 2013 as DOI: 10.4187/respcare.02786
AARC	Clinical Practice Guideline „Capillary Blood Gas Sampling for Neonatal & Pediatric Patients" - Respiratory Care (Respir Care 2001;46(5):506–513)
AARC	Clinical Practice Guideline „Sampling for Arterial Blood Gas Analysis" - Respiratory Care (Respir Care 1992;37:913-917)

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Literatür	
IFCC	Approved IFCC Recommendations on Whole Blood Sampling, Transport and Storage for Simultaneous Determination of pH, Blood Gases and Electrolytes" Eur J Clin Chem Clin Biochem. 1995 Apr;33(4):247-53.
AARC	Clinical Practice Guideline "Blood Gas Analysis and Hemoximetry: 2013", published on July 30, 2013 as DOI: 10.4187/respcare.02786
AARC	Clinical Practice Guideline "Capillary Blood Gas Sampling for Neonatal & Pediatric Patients" - Respiratory Care (Respir Care 2001;46(5):506–513)
AARC	Clinical Practice Guideline "Sampling for Arterial Blood Gas Analysis" - Respiratory Care (Respir Care 1992;37:913-917)

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